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
How to Purchase:
Contact Becker Orthopedic
(800) 521-2192 toll-free
(248) 588-7480

For a list of distributors worldwide:
www.beckerorthopedic.com

©2023 Tamarack Flexure Joint is a registered trademark of Tamarack Habilitation Technologies, Inc.
Made in the USA
REV 8/2023

Tamarack® Plantarflexion Motion Limiter Kit



Medial and Lateral,
Mold-In Adjustable Stops (1 Pair)

Part #741-ML-PF

www.tamarackhti.com



MADE IN MINNESOTA ★ U.S.A. ★

Fabrication & Installation Instructions

- To ensure that the adjustable stop hardware will incorporate into the AFO shell during molding, stockinette must not be placed over the motion limiter (joint stop) components.
Using brads provided, anchor the plantarflexion stop hardware posterior of the Tamarack Flexure Joint molding dummy. The posterior distance should be a minimum of 1" (25.4mm) from joint center line, which is the length of the molding spacer tab provided. A greater distance between the molding dummy and joint stop will work fine, but result in a smaller range of adjustability.
- Vacuum-form per standard procedure.
- Once plastic has sufficiently cooled and removed from mold, remove Flexure Joint Dummy from AFO; leaving motion limiter stop in place (see photo II).
- Separate AFO through centerline of Flexure Joint cavity and motion limiter using FINE TOOTH COPIGN SAW BLADE. (see photo III).
- Mark center of joint cavity to create reference point and ensure you don't grind here (see photo IV). Mark both foot and calf sections.
- Grind carefully to expose adjustment screw heads (see photo V). Use caution to avoid removing centerline mark on Flexure Joint Cavity.
- Grind to black stops, as needed for desired ankle range of motion, between 0° and 15° (see photo VI).
- Remove adjustment screw from adjustment screw housing (2).
- Use 5mm drill bit to pass through adjustment screw housing and drill hole on proximal side of motion limiter cavity (See photo VII)
- Re-insert set screw with head proximal. (See Photo VIII)
- Finish assembly and install Tamarack Flexure Joints (see photo IX).
- When adjustment screws are in their final position, per patient gait evaluation, use thread locker to hold in place.

Photo below depicts finished product

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