



Rivet Cover Patch


*Also Available in precut patch
kits from your O&P supplier*




Oval Patch
MED & LRG




Metatarsal Patch (right)
SML, MED & LRG




Circle Patch
SML, MED & LRG



Metatarsal Patch (left)
SML, MED & LRG



Leaf Patch
SML, MED & LRG



Teardrop Patch
SML, MED & LRG




Patch Application Recommendations



Bar Patch

Ideal shape for above & below-knee prosthetic socket and AFO brims



Leaf Patch

Ideal for prosthetic socket walls and concave surfaces of foot orthoses. Notches around perimeter prevent unwanted wrinkles on patch surface



Oval Patch

Multi-purpose shape for prosthetic sockets, AFOs, footwear, insoles, and spinal orthoses



Arch (Heel) Patch

Ideal for patching worn-out shoe linings on athletic and therapeutic footwear

For additional ShearBan application tips & techniques,
download the ShearBan Product Guide Online at www.shearban.com