

Toe Patch (left)

SML, MED & LRG

Toe Patch (right)

SML, MED & LRG

Bar Patch

MED & LRG

Arch Patch

SML, MED & LRG


Rivet Cover Patch

*Also Available in precut patch
kits from your O&P supplier*




SCAN TO VIEW THIS DOCUMENT ONLINE


as well as application tips, techniques,
& instruction guides for ShearBan® and
other Tamarack® products!




Oval Patch
MED & LRG




Metatarsal Patch (right)
SML, MED & LRG




Circle Patch
SML, MED & LRG



Metatarsal Patch (left)
SML, MED & LRG



Leaf Patch
SML, MED & LRG




Teardrop Patch
SML, MED & LRG

Patch Application Recommendations



Bar Patch

Ideal shape for above & below-knee prosthetic socket and AFO brims



Leaf Patch

Ideal for prosthetic socket walls and concave surfaces of foot orthoses. Notches around perimeter prevent unwanted wrinkles on patch surface



Oval Patch

Multi-purpose shape for prosthetic sockets, AFOs, footwear, insoles, and spinal orthoses



Arch (Heel) Patch

Ideal for patching worn-out shoe linings on athletic and therapeutic footwear



SCAN TO VIEW THIS DOCUMENT ONLINE

as well as application tips, techniques, & instruction guides for ShearBan[®] and other Tamarack[®] products!