




Tamarack[®]
HABILITATION TECHNOLOGIES, INC.

MADE IN MINNESOTA ★ U.S.A. ★

GLIDE WEAR[®] by Tamarack

Introducing an Innovative Solution to Skin Trauma for Wheelchair Users

Dual-layer Shear Reduction Zone™
glides with users skin.



Single layer stability zone enhances
user positioning and control

The GlideWear[®] Seating Interface
slips over most wheelchair cushion models,
providing protection from friction & shear.

For more information, visit www.GlideWear.com



APPLICATIONS

GlideWear® is indicated for wheelchair users who are at risk for skin breakdown due to friction and shear. GlideWear® technology increases the margin of safety against skin breakdown and complements skin protection & positioning cushion technologies.

SPECIFICATIONS

Construction: 85% Nylon, 15% Spandex. 100% Polyester Thread. Meets 16CFR1610 standards.

Sizes: Choose from 3 sizes, based on current cushion size.

Small: Fits cushions size 14" X 14" to 16" X 16"

Medium: Fits cushions size 16" X 16" to 18" X 18"

Large: Fits cushions size 18" X 18" to 20" X 20"

Warranty: Performance concerns and manufacturing defects are warranted for 12 months from date-of-purchase.

FEATURES & BENEFITS

Shear Reduction Zone™

Dual-layer fabric technology designed to slide against itself, absorbing friction & shear forces under bony prominences to assure a healthy margin of safety where skin is most susceptible to injury.

Positioning & Stabilization

Friction is retained in the areas surrounding the Shear Reduction Zone™, with a unique design within the interface, which helps maintain proper positioning and stability while sitting.

Versatile

GlideWear® slip-on cover is designed to fit most wheelchair cushions and may be installed over the existing cushion cover.

Breathable

Lightweight, breathable fabric helps reduce moisture buildup to preserve skin integrity

**Contact your seating & mobility provider to determine whether
GlideWear® is appropriate for your seating & skin health needs.**

**GlideWear® Technology is also featured on Comfort Company
and Vicair® brand wheelchair cushions.**

GlideWear® Technology Developed by



Tamarack®
HABILITATION TECHNOLOGIES, INC.